Lampade a fessura SL 980 Slit lamps

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
*The modern optic project included among the characteristics of the **SL 980** series is based on a system of multilayered antireflection treatment that ensures a more effective transmission of light thus giving both optic resolution and contrast an improvement of 20%. The optic quality produces a 20% increase in optic transmission and a 10% increase in depth of field.*

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*The **SL 980 Zoom** series has a zoom stereoscopic microscope with 8° convergent curve of the internal tubes, equipped with standard extractable eyepieces. Continuous enlargements from 6x to 33x, giving clear, brilliant and well contrasted images thanks to the multilayered antireflection treatment.*

On reference enlargement (red mark on lever for changing enlargement) can be used micrometer eyepieces.

Can be also used all models of beam splitter and accessories suited for the Galileian microscope.

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
The Greenough stereoscopic microscope 2x, thanks to the 13° convergence of binocular tube, and the Galileian stereoscopic microscope 3x - 5x and Zoom, thanks to the 8° convergence of binocular tube, combined to the Galileian system, allow an improved image fusion reducing the condition of stress for the observer.

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Illumination

The 6V 20W halogen bulb assures high quality color temperature. Available upon request the tilting device (prisma-holder head).

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Yellow filters insertion for fluoresceinic test

The CSO slit lamps on Galileian microscopes with 3 and 5 enlargements and Zoom microscopes allow the option of having the yellow filters OG 530 (fluoresceinic test) directly inserted into the microscope box for faster testing and improved optic yield of the set.

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Light source

Halogen illumination is transmitted more effectively by means of a new optic system that produces a clear and brilliant image. A 17% increase in aperture (Ø 14 mm. on all models).

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

A single joystick controls elevation and horizontal position. A rheostat for the luminous intensity adjustment of the slit lamp (optional) is fitted near the joystick. As a result operations are simple, fast and extremely precise.

Microscopio Galileiano con microcamera montata
Galileian microscope with microcamera

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
Beam splitter insertion

Loosen the knob, separate the microscope unit and insert the beam splitter. Then remount the microscope unit on the beam splitter and tighten the knob.

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*The SL 980 5x
slit lamp mounting the Z 800
applanation tonometer*

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Beam splitter

*system for optic splitting
of the microscope image
(50% on videocamera
- 50% on operator's eye).*

*It can be supplied in
the following models:*

- single plug for Panasonic™
microcameras
- single plug for C gauge
telecameras
- single plug (requires adaptor
joint for telecameras)

*- double plug
they are also easy to install
on Galileian and Zoom
microscopes.*

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Windows™ software - PC - digital monitor and colour printer.

The system permits images and digital video acquisition and elaboration.

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
Eye Image® is

- image and digital video acquisition, analysis and filing
- image measuring, editing and image printing
- patient examination files
- statistics
- personalized print and report
- unit expansion for corneal topography and applications of contact lenses.

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n° 96108701.2 ITB0950577

The system, with SL 980 Speedy, permits video shots without interfering with the observation optics, ensuring that operator and camera workings contemporary with 100% luminosity efficiency.


The camera infact uses for the video shots the objective that is not used by the observer.

This CSO technical solution is patented
n° 96108701.2 ITB0950577

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Optional accessories SL 980

- 16x eyepieces for 980 2x
- 10x eyepieces for 980 3x/5x
- Hruby lens kit
- Tonometer mount
- Hand-rest
- Rheostat on instrument base
- Video adapter
- Micrometer eyepiece
- Second observer tube
- Bearing with contrast and polarization filters
- Beam splitter for microcamera
- Beam splitter with one plug only for standard videocamera
- Beam splitter with two plugs (reflex camera-second observer-video)
- Light diffusion on prisma-holder head
- Yellow filters included in microscope
- Joystick's push button (useful in Eye Image® acquisition).

SL 980 standard accessories


- Chin-rest tissue
- Spare illumination bulb
- Fuses
- Maintenance kit

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Yellow or
polarising filters
device

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The diagnostic functionalities of the SL 980 slit lamps are constantly developed thanks to the collaboration from surgical and microsurgical divisions of the greatest hospitals, universities and research centres, medical institutions, clinics and eye centres.

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© C.S.O. srl 2004. I marchi citati appartengono ai rispettivi proprietari - All trademarks mentioned are of their respective companies
C.S.O. S.r.l. Costruzione Strumenti Oftalmici via degli Stagnacci,12/E